

eBite

Wireless & Rechargeable
Intraoral Lighting System


Liberating dentistry, one hand at a time



Wireless



Rechargeable




Bite Block



Light



■ Product Features



■ Specifications

Biteblock	Autoclavable Silicone material (Large and small sizes available)
LED Light	Illuminance: 5,000~8,000 lux Color temperature: 5,500K Operation life: 50,000 hours
Battery	Li-Polymer Rechargeable battery
Power Adapter	Input: 100~200VAC / 50~60Hz Output: DC 5.0V / 1.0A