

# eBiteplus


Intraoral Suction & Lighting System

Liberating dentistry, one hand at a time



Suction

+



Bite Block

+



Tongue Retractor

+



Light



Free up your assistants' hands while securing valuable workspace during routine dental procedures.

The eBite Plus combines traditional tools to make procedures more efficient for dentists and more comfortable for their patients.



- Autoclavable hand piece, bite blocks, and tongue retractors
- Multiple sizes for bite blocks and tongue retractors
- On-off light switch with three different light intensity levels
- Compatible with both HVE and SE

■ Product Features



■ Specifications

Handpiece	Autoclavable
Biteblock & Tongue Retractor	Autoclavable Silicone material (Large and small sizes available)
LED Light	Illuminance: 5,000~8,000 lux Color temperature: 5,500K Operation life: 50,000 hours
Suction Pipe	Independent control system for Left and Right Stainless steel twin type
Suction Hose	OD: 0.472 inches ID: 0.236 inches
Power Adapter	Input: 100~200VAC / 50~60Hz Output: DC 4.5V / 2.0A